

COAXIAL MICROPHACO SET

Your step to the 2.0 mm Incision


6-20/6-107 Coaxial Microphaco Keratome Diamond Knife

- Angled Trapezoid Self-Diving Blade 2.0/2.3mm
- ThermoCut® Technology




6-10/6-0501 Micro SidePort Diamond Knife

- Straight 0.6mm 45° Diamond blade
- Ideal for making the SidePort during the Coaxial MicroPhaco procedure and flexible iris retractors implantation
- ThermoCut® Technology



4-0391S Inamura Capsulorrhexis Forceps

- Cystotome tips, 10mm jaws
- Tips do not open more than 2mm



7-1162S Akahoshi Combo Prechopper for Coaxial Microphaco

- Specially designed to split hard and soft nuclei through 2.0mm incision



16-2804 Tsar II Injector

- Enables a safe single-hand IOL implantation using the Counter Traction technique and C cartridge through 2.0mm incision



8601 4th Street North, Suite 201,
St. Petersburg, FL 33702, USA

Tel.: 727-568-0909

Toll Free: 877-77-RUMEX

Fax: 727-568-0919

E-mail: rumex@rumex.net



Shop @ www.rumex.net