



## Toric IOL Marker




**3-181**

Allows easy and precise marking of the steep meridian to position the Toric IOL. Specially designed pins enable easy rotation.




## LRI Marker



**3-180** 6 & 8 mm Chord Lengths

**3-1801** 40-60-80 Degrees

Intra-Op Marker  
Designed for Limbal Relaxing Incision




13770 58th Street North, Suite 317  
Clearwater, FL 33760, USA  
Tel.: 727-568-0909  
Toll Free: 877-77-RUMEX  
Fax: 727-568-0919

Hotline for customers in  
Europe, Asia,  
Africa & Latin America:  
Tel: +371 6616 3182  
e-mail: rumex@rumex.net

## Whitehouse Gravity Toric/LRI Marker *with hinder position of the weight*



- 3-193** *Angled*
- 3-1931** *Straight*
- 3-1932** *Straight (with reduced diameters)*



Pre-Op reference marker designed by *Dr. Geoffrey Whitehouse* to accurately mark the limbus at 3, 6 and 9 o'clock to avoid the effect of cyclotriorion.

## Toric IOL Combo Marker



- 3-194** *Vertical 0-0 axis*
- 3-1941** *Horizontal 0-0 axis*



Pre-Op marker to mark the axes while the patient is in sitting position at the slit lamp. Reduces cyclotorsion effect in supine position.